

Enteral Nutrition



Enteral Syringe

- FDA cleared for direct oral administration and enteral delivery
- The ENFit syringe tip is compliant with ISO 80369-3 standards
- ENFit syringes are often used in pharmacies, NICU/ PICU, adult ICU, alternate care, home care, and physician offices.

ENFit Single Use Syringe



Product Feature :

- ENFit female.
- DEHP Free
- Latex Free
- Low dose
- Amber color is available



ENFit Home Use O-Ring Syringe

- Up to 84 uses
- O-Ring Plunger
- Low resistance
- Amber color is available



Sizes	Low Dose	Single Use	Home Use
1ml	•	ES1	ESH1
3ml	•	ES3	ESH3
5ml	•	ES5	ESH5
10ml		ES10	ESH10
20ml		ES20	ESH20
35ml		ES35	ESH35
60ml		ES60	ESH60
100ml		ES100	N/A





Thumb Control Ring Syringe

- Designed for maximum control and one-handed administration
- Re-sealable pole bag designed for safe and convenient storage

Flat Top Syringe


- Visible dosing marks designed for ease of use

Sizes	Enfit Tip	Thumb Control Ring	Flat Top	Code
60ml		•		URTS60
60ml			•	URFS60
60ml	•	•		URES60
60ml	•		•	ES60
100ml	•		•	ES100
100ml		•		URTS100



The Enfit Connector

- The ENFit connector is compliant with the ISO80369-3.
- Supply you transition adaptor from Legacy connector to Enfit. (ENC3)
- ENFit syringe cap used as a protective cap to temporarily secure any remaining nutrition between feedings. (ENC6)
- Supply you transition adaptor from Enfit to Legacy connector. (ENC10)



ENFit Feeding Set

Feature:

- ENFit feeding sets are compliant with the ISO 80369-3 to help reduce the risk of misconnections.
- Latex-free and DEHP-free.
- EVA material available.
- G-Gravity use, P-pump use, ENF01-05 available G type and P type.



ENF05G
Universal Set



ENF01G
ENPLUS spike with cap



ENF02P
Feeding bag
Sizes: 500/1000/1200/1500ml



ENF03G
Feed cap adaptor



ENF04G
Spike with feeding connector




Oral Dispenser

Feature:

- Special tip design for wrong route administration prevention.
- O-Ring plunger design is the preferred option for smooth and accurate delivery.
- Amber barrel design to protect light-sensitive drug.




O-Ring Plunger



Type 1

Sizes	Type 1	Type 2
1ml	ORS1-01	ORS2-01
3ml	ORS1-03	ORS2-03
5ml	ORS1-05	ORS2-05
10ml	ORS1-10	ORS2-10
20ml	ORS1-20	ORS2-20

Type 2



ENFit Stepped

We're on the path to safer connections

The ENFit connector reduces the risk of tube feeding misconnections, helping enhance patient safety. It has a unique enteral-specific design.



ENC10

- Has a locking feature that secures the safety connection with administration sets and syringes that have ENFit connectors.
- Aligns with global, industry-wide standards specified by ISO and recommended by the FDA.
- Supply you transition adaptor from Enfit to Legacy connector